

# A

TRUREPAIR Solutions.....A-1

## TRUREPAIR<sup>®</sup> Solutions

### Product Highlights



#### TRUREPAIR Solutions

Smith & Nephew TRUREPAIR devices provide a resorbable scaffold and related instrumentation for the repair and replacement of bone. TRUREPAIR devices are made from proprietary POLYGRAFT<sup>®</sup> technology—a completely synthetic composite material that is biologically friendly and biomechanically stable.

The TRUREPAIR product family includes TRUFIT<sup>®</sup> BGS Plugs, TRUGRAFT<sup>®</sup> Syringes and Granules, TRUWEDGE<sup>®</sup> implants and TRUBLOCK<sup>®</sup> implants.

## Tissue Repair

### TRUREPAIR® Product line

The Smith & Nephew TRUREPAIR product line provides porous structural, bioresorbable scaffolds uniquely designed for the repair and replacement of bone.

Smith & Nephew is focused on developing Innovative Tissue Repair Solutions for sports medicine indications. Utilizing a proprietary combination of polymers, ceramics, and fibers, the material is designed to be a highly porous scaffold to support bone incorporation and remodeling, it's biomechanically sound during the time required for healing and it's biologically friendly by absorbing biologic fluids and nutrients. As innovators in tissue repair, Smith & Nephew will continue to develop products designed to reduce patient morbidity and improve patient outcomes.

Our products include the TRUFIT® BGS Plugs for use in filling bony voids or gaps caused by trauma or surgery that are not intrinsic to the stability of the bony structure. TRUKOR® Plus Instruments for the preparation of precise circular sites for the TRUFIT BGS Plug. TRUGRAFT® BGS Syringe to offer a more convenient means of delivering TRUGRAFT BGS Granules into bone void sites with limited accessibility.

### TRUGRAFT BGS Granules and Syringe



REF	Description
72200919	TRUGRAFT BGS Granules, 10 cc (2 g)
72200920	TRUGRAFT BGS Granules, 30 cc (6 g)
72200921	TRUGRAFT BGS Syringe, 5 cc (1 syringe per package)
72200922	TRUGRAFT BGS Syringe, 10 cc (2 syringe per package)

### TRUFIT BGS Plug Kits (Plug, Delivery Device, Trimming Knife)

The safety and effectiveness of this device for use in osteochondral defects have not been established.



REF	Description
72200923	TRUFIT BGS Plug – Size Code Green, 5 mm
72200924	TRUFIT BGS Plug – Size Code Yellow, 6 mm
72200925	TRUFIT BGS Plug – Size Code Red, 7 mm
72200926	TRUFIT BGS Plug – Size Code Tan, 7.5 mm
72200927	TRUFIT BGS Plug – Size Code Blue, 9 mm
72200928	TRUFIT BGS Plug – Size Code Orange, 9.5 mm
72200929	TRUFIT BGS Plug – Size Code Purple, 11 mm
72200948	TRUFIT Ruler

## TRUREPAIR<sup>◇</sup> Solutions

### TRUKOR<sup>®</sup> Plus Instruments



REF	Description
72200965	TRUKOR Plus Site Preparation Kit (Contains Sterilization Tray, Drills, Universal Drill Handle, Caps, Sizers, Depth Gauge, Obturators)
72200964	TRUKOR Plus Instrument Sterilization Tray
72200954	TRUKOR Plus Drill – Size Code Green, 5 mm
72200955	TRUKOR Plus Drill – Size Code Red, 7 mm
72200956	TRUKOR Plus Drill – Size Code Blue, 9 mm
72200957	TRUKOR Plus Drill – Size Code Purple, 11 mm
72200962	TRUKOR Plus Universal Drill Handle
72200941	TRUKOR Drill Sleeve Cap – Size Code Green, 5 mm
72200942	TRUKOR Drill Sleeve Cap – Size Code Red, 7 mm
72200943	TRUKOR Drill Sleeve Cap – Size Code Blue, 9 mm
72200944	TRUKOR Drill Sleeve Cap – Size Code Purple, 11 mm
72200949	TRUKOR Obturator – Size Code Green, 5 mm
72200950	TRUKOR Obturator – Size Code Red, 7 mm
72200951	TRUKOR Obturator – Size Code Blue, 9 mm
72200952	TRUKOR Obturator – Size Code Purple, 11 mm
72200958	TRUKOR Plus Sizer – Size Code Green, 5 mm
72200959	TRUKOR Plus Sizer – Size Code Red, 7 mm
72200960	TRUKOR Plus Sizer – Size Code Blue, 9 mm
72200961	TRUKOR Plus Sizer – Size Code Purple, 11 mm
72200963	TRUKOR Plus Depth Gauge

### TRUKOR<sup>®</sup> Instruments - Drill Sleeves



REF	Description
72200937	TRUKOR Drill Sleeve – Size Code Green, 5 mm, disposable
72200938	TRUKOR Drill Sleeve – Size Code Red, 7 mm, disposable
72200939	TRUKOR Drill Sleeve – Size Code Blue, 9 mm, disposable
72200940	TRUKOR Drill Sleeve – Size Code Purple, 11 mm, disposable

## TRUREPAIR<sup>◇</sup> Solutions

### TRUBLOCK<sup>◇</sup> BGS



REF	Description
72200967	TRUBLOCK BGS – 15 mm X 15 mm X 20 mm
72200968	TRUBLOCK BGS – 20 mm X 20 mm X 20 mm
72200969	TRUBLOCK BGS – 25 mm X 25 mm X 20 mm

### TRUWEDGE<sup>◇</sup> BGS (Length 40 mm, Width 15 mm)



72200970	TRUWEDGE BGS – 7.5 mm (Height), 7.5°
72200971	TRUWEDGE BGS – 10 mm (Height), 10°
72200972	TRUWEDGE BGS – 12.5 mm (Height), 12.5°
72200973	TRUWEDGE BGS – 15 mm (Height), 15°
72201026	TRUWEDGE BGS – 5 mm (Height), 5°

The TRUFIT<sup>◇</sup> CB products are indicated for cartilage and bone repair. These products are *not* available for distribution in the United States and its territories.

### TRUFIT CB Plug Kits (Plug, Delivery Device, Trimming Knife)



72200930	TRUFIT CB Plug – Size Code Green, 5 mm
72200931	TRUFIT CB Plug – Size Code Yellow, 6 mm
72200932	TRUFIT CB Plug – Size Code Red, 7 mm
72200933	TRUFIT CB Plug – Size Code Tan, 7.5 mm
72200934	TRUFIT CB Plug – Size Code Blue, 9 mm
72200935	TRUFIT CB Plug – Size Code Orange, 9.5 mm
72200936	TRUFIT CB Plug – Size Code Purple, 11 mm

### TRUFIT CB Angle Plugs (Plug, Delivery Device, Trimming Knife)



72201014	TRUFIT CB Angle Plug, Size Code Green, 5 mm
72201015	TRUFIT CB Angle Plug, Size Code Red, 7 mm
72201016	TRUFIT CB Angle Plug, Size Code Blue, 9 mm

### TRUFIT CB Saddle Plugs (Plug, Delivery Device, Trimming Knife)

72201909	TRUFIT CB Saddle Plug, Size Code Blue, 9 mm
72201910	TRUFIT CB Saddle Plug, Size Code Purple, 11 mm

