

B

Arthroscopes

VideoArthroscopes, Autoclavable	B-1
VideoArthroscopes, Non-Autoclavable	B-2
Direct-View Arthroscopes	B-4
Direct-View Arthroscopes, Autoclavable	B-4
Direct-View Arthroscopes, Non-Autoclavable	B-4
Specialty Endoscopes (ECTRA® II) and Accessories	B-5
Made to Order	B-6

VideoArthroscopes, Autoclavable



Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Wide-angle lens design for a larger field-of-view

Unparalleled brightness

Sharp images throughout the entire field-of-view

Lightweight construction

Sapphire lenses for improved scratch resistance over conventional glass

Epoxy-free proprietary seal design guarantees leak-free multiple sterilization modality

Compatible with Smith & Nephew autoclavable Camera Heads

Validated for both steam and chemical (EtO, Steris®, and Sterrad®) sterilization modalities

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
<i>VideoArthroscopes, Autoclavable, High Definition 560H only</i>							
72200312	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
72200089	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
<i>VideoArthroscopes, Autoclavable, 470H head, Smith & Nephew mount</i>							
7210716	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
7209401	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm
7209402	2.7 mm	30°	90°	120 mm	Opposite	J-Lock	4.0 mm
7205929	2.7 mm	30°	90°	67 mm	Opposite	J-Lock	4.0 mm

REF	Description
7205548	Aluminum oxide vial for cleaning autoclavable endoscope lenses
7210721	Protective cap for autoclavable scopes. Fits 4.0 mm and 5.5 mm scopes. Package of 70.

Note: These VideoArthroscopes for use only with Smith & Nephew Autoclavable Camera Heads

VideoArthroscopes, Non-Autoclavable



Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Stainless steel construction to withstand chemical sterilization modalities (EtO, Steris®, and Sterrad®)

Wide angle of view

Ergonomically designed to directly attach to the camera head for ease of use

Insertion needle is fully rotatable without the need to reposition the camera head

Compatible with Smith & Nephew Non-Autoclavable Camera Heads (C-Mount), as listed below

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
-----	----------------	-------------------	---------------	----------------	----------------------	--------------------	-------------------------

Black Focus Ring, 19.5 mm focal length

Note: These VideoArthroscopes, non-autoclavable HD for use only with the Smith & Nephew 560H camera heads only.

72201501	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
72201671	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm
72201689	4.0 mm	30°	115°	160 mm	Same Side	J-Lock	5.0 mm

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
-----	----------------	-------------------	---------------	----------------	----------------------	--------------------	-------------------------

Blue Focus Ring, 35 mm focal length

7205922*	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
3766	4.0 mm	30°	115°	160 mm	Same Side	J-Lock	5.0 mm
3628	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm
3632	2.7 mm	30°	90°	120 mm	Opposite	J-Lock	4.0 mm
3768	2.7 mm	30°	90°	67 mm	Opposite	J-Lock	4.0 mm
3863	1.9 mm	30°	75°	65 mm	Opposite	J-Lock	2.0 mm

Note: These VideoArthroscopes for use only with the Smith & Nephew ED3, 360, 460, HD900, and 111 camera heads.
Not compatible with the Smith & Nephew 470 or 560 camera heads

REF	Description
8001701	Replacement O-ring for C-mount Video Arthroscopes

VideoArthroscopes, Non-Autoclavable



Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx
 Stainless steel construction to withstand sterilization modalities (EtO, Steris®, and Sterrad®)
 Wide angle of view
 Ergonomically designed to directly attach to the camera head for ease of use
 Insertion needle is fully rotatable without the need to reposition the camera head
 Compatible with Smith & Nephew Non-Autoclavable Camera Heads (C-Mount)

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
4125	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
4126	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm
4130*	2.7 mm	30°	90°	120 mm	Opposite	J-Lock	5.0 mm
4131†	2.7 mm	30°	90°	67 mm	Opposite	J-Lock	5.0 mm
4132†	2.7 mm	70°	90°	70 mm	Opposite	J-Lock	5.0 mm
7208133†	1.9 mm	30°	75°	62 mm	Opposite	J-Lock	2.0 mm


REF	Description
8001701	O ring

* Item shown

† HD and non-HD compatible

Direct-View Arthroscopes

Direct-View Arthroscopes, Autoclavable



Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Fully autoclavable, with both distal and proximal ends soldered for longevity

Wide angle of view

Autoclavable arthroscopes have sapphire distal lenses and rear ocular windows for much greater scratch resistance than conventional glass

Stainless steel construction to withstand sterilization modalities (EtO, Steris®, Sterrad®, and autoclave)

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
72202087	4.0 mm	30°	115°	160 mm	Opposite	J-Lock	5.0 mm
72202088	4.0 mm	70°	115°	160 mm	Opposite	J-Lock	5.0 mm
72202089	4.0 mm	45°	115°	157 mm	Opposite	J-Lock	5.0 mm
72202090	4.0 mm	0°	115°	160 mm	Opposite	J-Lock	5.0 mm
7208132	4.0 mm	30°	115°	175 mm	Opposite	Storz	5.0 mm
72200733	2.7 mm	70°	80°	119 mm	Opposite	J-Lock	5.0 mm
72201747	2.7 mm	30°	80°	67 mm	Same side	J-Lock	5.0 mm
7205682	2.7 mm	30°	90°	67 mm	Opposite	J-Lock	5.0 mm
7205681	2.7 mm	30°	80°	120 mm	Opposite	J-Lock	5.0 mm

Note: Compatible with all camera head with couplers

Direct-View Arthroscopes, Non-Autoclavable

Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Wide angle of view

Stainless steel construction to withstand sterilization modalities (EtO, Steris®, and Sterrad®)

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
4184	1.9 mm	30°	80°	67 mm	Opposite	J-Lock	2.0 mm

REF	Description
-----	-------------

7205548	Aluminum oxide vial for cleaning autoclavable endoscope lenses
7210721	Protective cap for autoclavable scopes. Fits 4.0 mm and 5.5 mm scopes. Package of 70.

Specialty Endoscopes (ECTRA[◊] II) and Accessories



Indications: Endoscopic release of carpal ligament, plantar fascia, and the extra-articular lateral ligament

Wide angle of view

Stainless steel construction to withstand sterilization modalities (EtO, Steris[®], and Sterrad[®])

REF	Type	Focal Length	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Recommended Light Cable
4144*	Direct-View	N/A	4.0 mm	30°	100°	106 mm	Same Side	5.0 mm
3896	Video	35 mm	4.0 mm	30°	100°	106 mm	Same Side	5.0 mm
72201597	Direct-View	N/A	4.0 mm	30°	115°	278 mm	Opposite Side	5.0 mm
72201598	Direct-View	N/A	4.0 mm	70°	115°	280 mm	Opposite Side	5.0 mm

Replacement Parts

Light post Adaptor Sleeve: Fits all Arthroscopes, Hysteroscopes/Cystourethroscopes, and Laparoscopes with 4 mm light post



REF	Description
7207900	For Storz [®] Light Guide, outside threads, full length



4347*	For Smith & Nephew/Wolf [®] Light Guide
-------	--

8001701	Replacement O-ring for videoscope and C-mount couplers
---------	--

* Item shown

Made to Order

Made to Order products have longer production times. Please contact your sales representative or customer service for more information.

VideoArthroscope, Multimode

Compatible with Smith & Nephew Non-Autoclavable Camera Heads (C-Mount)

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
7209911Z	4.0 mm	30°	115°	160 mm	Same Side	J-Lock	5.0 mm
7210750	4.0 mm	30°	115°	97 mm	Same Side	J-Lock	5.0 mm

VideoArthroscopes, non-autoclavable

Black Focus Ring, 25 mm focal length

Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Stainless steel construction to withstand sterilization modalities (EtO, Steris®, and Sterrad®)

Wide angle of view

Ergonomically designed to directly attach to the camera head for ease of use

Insertion needle is fully rotatable without the need to reposition the camera head

Compatible with Smith & Nephew Non-Autoclavable Camera Heads (C-Mount)

Not for use with HD camera heads

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
4486	4.0 mm	30°	115°	160 mm	Same Side	J-Lock	5.0 mm

VideoArthroscopes, Non-Autoclavable

Blue Focus Ring, 35 mm focal length

Indications: Arthroscopy and examination/treatment of nasal cavities and nasal pharynx

Stainless steel construction to withstand chemical sterilization modalities (EtO, Steris®, and Sterrad®)

Wide angle of view

Ergonomically designed to directly attach to the camera head for ease of use

Insertion needle is fully rotatable without the need to reposition the camera head

Compatible with Smith & Nephew Non-Autoclavable Camera Heads (C-Mount)

Not for use with HD camera heads

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
7210987	1.9 mm	30°	75°	65 mm	Opposite	J-Lock	5.0 mm
3797	2.7 mm	70°	90°	70 mm	Opposite	J-Lock	5.0 mm

Made to Order

Made to Order products have longer production times. Please contact your sales representative or customer service for more information.

Direct-View Arthroscopes, Autoclavable

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
7210694	2.7 mm	0°	80°	67 mm	Opposite	J-Lock	5.0 mm
7210910	2.7 mm	0°	80°	120 mm	Opposite	J-Lock	5.0 mm
7204734	4.0 mm	0°	115°	175 mm	Opposite	Storz°	5.0 mm
7204735	4.0 mm	70°	115°	175 mm	Opposite	Storz	5.0 mm

Direct-View Arthroscopes, Non-Autoclavable

REF	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Cannula Connection	Recommended Light Cable
3948	2.7 mm	70°	90°	70 mm	Opposite	J-Lock	5.0 mm

Specialty Endoscopes

REF	Type	Focal Length	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Recommended Light Cable
4139Z	Video	25 mm	4.0 mm	30°	115°	106 mm	Same Side	5.0 mm

Diagnostic Cannulas, Extended Length

For use with hip arthroscopy scopes REF 7211054 and 7211055

REF	Description
7211075	5.8 mm diagnostic cannula, single valve, rotatable
7207913	5.8 mm diagnostic cannula, double valve, rotatable
7207898	Cannulated obturator for cannulas 7211075 and 7207913

