

D

Smith & Nephew DYONICS® POWER
Control SystemsD-1
Smith & Nephew DYONICS® Motor
Drive UnitsD-2
AccessoriesD-3
Replacement PartsD-4

Shaver Systems

Product Highlights



DYONICS POWERMINI
Small Joint Handpiece
Page D-2

NEW! Smith & Nephew DYONICS POWERMINI Small Joint Handpiece with Blade Multi-Positioning, hand controls, suction lever.



DYONICS POWER II Control Unit
Page D-2

DYONICS POWER II Control Unit
Combines the most aggressive arthroscopic soft tissue resection with powerful torque for efficient bone removal. Compatible with all DYONICS POWER handpieces and footswitches.



DYONICS POWER II Footswitch
Page D-3

DYONICS POWER II Footswitch
Ergonomic design may be operated in variable speed or on/off modes. Provides flexibility to operate in variable speed or on/off modes and activate handpiece and footswitch controls during a procedure.

Smith & Nephew DYONICS[®] POWER Control Systems



Smith & Nephew DYONICS POWER Shaver Systems

The following features are standard on all DYONICS POWER control units:

- Blade family auto-recognition recalls blade and burr settings
- Blade window lock helps protect tissue from open edgeforms
- Communicates with DYONICS 25 and InteliJET[®] pumps to maximize visibility while minimizing fluid usage.

In addition, the high-performance DYONICS POWER II Control Unit offers two unique, customizable oscillate modes for efficient soft tissue resection; dual hot ports for simultaneous shaver and powered instrument operation; and an integrated flush lavage feature to rapidly clear the visual field of debris.

REF	Description
72200873	Smith & Nephew DYONICS POWER II Shaver System Control Unit. All units 100-240 VAC, 50/60 Hz
72201030	U.S. system
72201034	U.K. system
72201037	German system
72201040	Spanish system
72201043	French system
72201046	Italian system
72201049	Swedish system
72201052	Dutch system
72201055	Portuguese system
72201058	Danish system
72201061	Norwegian system
72201322	Korean system

All kits include a control unit, power cord, and operations manual.

REF	Description
7205841	Smith & Nephew DYONICS POWER Shaver System Control Unit. All units 100-240 VAC, 50/60 Hz
7205886	U.S. kit
7205887	U.K. kit
7205888	International kit
7205889	French kit
7205890	German kit
7205891	Dutch kit
7205892	Italian kit
7205893	Spanish kit
7205894	Swedish kit

All kits include a power cord, operations manual and communications cable for the DYONICS 25 Pump

Call +1-800-343-5717 (U.S.) or contact your local representative.

Smith & Nephew DYONICS[◇] Motor Drive Units

For use with Smith & Nephew DYONICS POWER and EP-1[°] Control Units

REF	Description
72200616*	Smith & Nephew DYONICS POWERMAX [°] ELITE Handpiece with blade multi-positioning, hand controls, suction lever, 10 ft (3 m) cord
72200617*	Smith & Nephew DYONICS POWERMAX ELITE Handpiece with Blade Multi-Positioning, no hand controls, suction lever, 10 ft (3 m) cord
72200872	Smith & Nephew DYONICS POWERMAX ELITE Handpiece hand-controlled with opposite-side suction lever

*NOTE: MDU's accept all blades

72201500	POWERMINI Small Joint Handpiece with blade multi-positioning, hand controls, suction lever, 10 ft. (3m) cord For use with Smith & Nephew DYONICS POWER II Control System.
72201503	POWERMINI Small Joint Handpiece with Blade Multi-Positioning, no hand controls, suction lever, 10 ft. (3m) cord For use with Smith & Nephew DYONICS POWER II Control System.

**NOTE: MDU only works with the "Disposable Small Joint Blades/Burrs" listed on page E-7.

7205357	Smith & Nephew DYONICS Mini-Motor Drive, suction lever, 10 ft (3 m) cord
7205971***	Smith & Nephew DYONICS POWER Ultralight Motor Drive, suction lever, 10 ft (3 m) cord

***NOTE: MDU accepts all straight blades, but only ULTRA series curved blades.

Accessories

Smith & Nephew DYONICS® POWER Footswitches



Low profile, cutting mode and blade window lock control, 12 ft (3.6 m) cord

REF	Description
7205399*	Smith & Nephew DYONICS POWER Footswitch, low profile, on/off
7205397	Smith & Nephew DYONICS POWER Footswitch, low profile, variable speed



Pedal-style, 12 ft (3.6 m) cord

REF	Description
7205396	Smith & Nephew DYONICS POWER Footswitch, pedal-style, on/off



DYONICS POWER II Footswitch

REF	Description
72201092	Pedal-style, on/off or variable speed. Operates with DYONICS POWER II Control. (REF. 72200873 only)

Sterilization Tray

REF	Description
4323	Shaver Sterilization Tray for Smith & Nephew DYONICS POWER shaver handpieces
4324	Shaver Sterilization Insert for Smith & Nephew DYONICS POWER shaver handpieces

* Item shown

Replacement Parts

Replacement Motor Drive Cord



REF	Description
7208207	Replacement Motor Drive Cord for Smith & Nephew DYONICS® POWER and EP-1® handpieces, standard or hand control type

***NOTE:** For standard and POWERMAX MDU's use 8101720 O Ring.

Replacement power cords and connectors

REF	Description
7205081	AC Power Cord for Smith & Nephew DYONICS POWER and EP-1 Shaver Control Units, 3 ft (90 cm), 120 VAC
4428	AC Power Cord for Smith & Nephew DYONICS POWER and EP-1 Shaver Control Units, 12 ft (3.6 m), 120 VAC
4476	AC Power Cord for Smith & Nephew DYONICS POWER and EP-1 Shaver Control Units, 3.25 ft (1 m), 230 VAC
7205393	Handpiece Connector, Smith & Nephew DYONICS POWER and EP-1 Control Unit, replacement
7205394	Footswitch Connector, Smith & Nephew DYONICS POWER and EP-1 Control Unit, replacement

Replacement Caps and Lanyards

The table below will assist you in ordering the replacement cap and lanyards. There are three (3) shaver handpiece motor drive cords, therefore, it is important that you reference the serial number on the motor drive cord to the table below to ensure that the proper replacement cap and lanyard is ordered.

System	Description	Replacement REF	Cable Serial Numbers	IFU
S&N DYONICS® POWER and EP-1®	MDU motor cord (SCREW-ON TYPE)	7205902 8100001	TJ 5524 and HIGHER	1060809
S&N DYONICS POWER and EP-1	MDU motor cord (SCREW-ON TYPE)	7205903 8101765	TJ 2860 to TJ 5523	
S&N DYONICS POWER and EP-1	MDU motor cord (PUSH-ON CAP)	8101544*	TJ 0000 to TJ 2859	
S&N DYONICS POWER and EP-1	FOOTSWITCH cord	8100002	ALL	
PS3500EP®	MDU motor cord	8101724	ALL	
Ultralight	MDU motor cord (SCREW-ON TYPE)	7205902 8100001	TJ 5524 and HIGHER	1060809

Replacement Motor Drive Cords

7208208	Replacement Motor Drive Cord for Smith & Nephew DYONICS POWER and EP-1 Mini Motor Drive
7209071	Replacement Motor Drive Cord for Smith & Nephew DYONICS POWER Ultralight Motor Drive

*NOTE: Soak cap 8101544 does not include lanyard.