

# E

Smith & Nephew DYONICS<sup>®</sup>  
ELECTROBLADE<sup>®</sup> Resector .....E-1

Smith & Nephew DYONICS<sup>®</sup>  
BONECUTTER<sup>®</sup> Blades .....E-2

Smith & Nephew DYONICS<sup>®</sup>  
ORBIT<sup>®</sup> Curved Blades .....E-3

Pre-Curved and Long Disposable Blades ...E-4

Straight Disposable Blades .....E-5

Disposable Burrs .....E-7

Disposable Small Joint Blades/Burrs .....E-8

Disposable Mini Blades/Burrs .....E-9

Disposable Mini Blades —  
Long and Curved .....E-10

Short Sheath Components .....E-11

Smith & Nephew offers a complete line of high quality, single-use blades and burrs for all arthroscopic procedures. Smith & Nephew has a wide range of blade styles and components, including curved blades, ORBIT blades, long blades, burrs, small joint blades, and short sheath and obturator handles.

Use of reprocessed, single-use blades may permanently damage, impede performance, or cause failure of your Arthroscopic Shaver System. Use of such products may render any warranties null and void.

## Sterile Disposable Blades and Burrs

### Product Highlights



#### INCISOR<sup>®</sup> PLUS ELITE blades

Thanks to a precise sharpening process, INCISOR PLUS ELITE blades deliver twice the cutting efficiency of competitive aggressive blades.

INCISOR<sup>®</sup> PLUS ELITE blades  
Page E-4



#### DYONICS BONECUTTER<sup>®</sup> ELECTROBLADE<sup>®</sup> Resection

The first all-in-one device for subacromial decompression, simultaneously resects and coagulates soft tissue and cuts bone.

DYONICS BONECUTTER<sup>®</sup>  
ELECTROBLADE<sup>®</sup> Resection  
Page E-2

#### DYONICS ELITE Burrs

Combine power, speed and responsiveness with precision sharpness for aggressive bone resection.



DYONICS ELITE Burrs  
Page E-9

## Smith & Nephew DYONICS<sup>◇</sup> ELECTROBLADE<sup>◇</sup> Resector



DYONICS ELECTROBLADE Resectors deliver simultaneous mechanical resection and coagulation of soft tissue. This helps increase procedure efficiency and maintain visual clarity throughout the procedure. In addition, the bipolar design eliminates the need for a ground pad. The new BONECUTTER ELECTROBLADE Resector is the first all-in-one device for subacromial decompression. It resects and simultaneously coagulates tissue and also cuts bone.

### DYONICS ELECTROBLADE Resectors

Packaged 3 per box, sterile

REF	Description	Latch Color
7205962	DYONICS BONECUTTER <sup>◇</sup> ELECTROBLADE Resector <i>Compatible with the VULCAN<sup>◇</sup> Generator (requires the VULCAN Generator Adapter) and "the Erbe Generator"</i>	Orange
7209855	VULCAN Generator Adapter <i>Required to operate 7205962, packaged non-sterile, reusable</i>	

REF	Description	Latch Color
7205961	4.5 Full Radius, Valleylab Generator-compatible <i>Compatible with the Valleylab*, Erbe** and VULCAN (when used with the VULCAN Generator Adapter) Generators</i>	Yellow
7209855	4.5 Full Radius, VULCAN <sup>◇</sup> Generator-compatible	Yellow
7209700	4.5 ELITE <sup>◇</sup> , Full Radius, Valleylab Generator-compatible <i>Compatible with the Valleylab*, Erbe** and VULCAN (when used with the VULCAN Generator Adapter) Generators</i>	Maroon
7209983	4.5 ELITE, VULCAN Generator-compatible <sup>†††</sup>	Maroon

\*Valleylab Generator models are the Valleylab Force FX(tm) and Force(tm) 2;

\*\*Erbe Generators are the Erbe 300 and 350

## Smith & Nephew DYONICS<sup>◇</sup> BONECUTTER<sup>◇</sup> Blades



### Smith & Nephew DYONICS BONECUTTER Blades


Packaged 6 per box, sterile

Resects both soft tissue and bone with a single blade, eliminating the need to switch from a straight blade to a burr

Proprietary SULILLIUM<sup>◇</sup> Coating helps maintain blade sharpness and durability throughout the procedure

Ideal for Notchplasty and Acromioplasty procedures

Operates at up to 8000 rpm for effective bone burring



REF	Description	Latch Color
7210486	3.5 Full Radius BONECUTTER Blade	Beige
7206011	4.5 Full Radius BONECUTTER Blade	Yellow
7206010	5.5 Full Radius BONECUTTER Blade	Orange
7206013	4.5 SYNOVATOR <sup>◇</sup> with a Full Radius inner BONECUTTER Blade	Forest green
7206012	5.5 SYNOVATOR with a Full Radius inner BONECUTTER Blade	Melon yellow

### Smith & Nephew DYONICS BONECUTTER Long Shaver Blades

18 cm working length, packaged 3 per box, sterile

REF	Description	Latch Color
7210979	4.5 Full Radius BONECUTTER Blade	Clear

## Smith & Nephew DYONICS<sup>®</sup> ORBIT<sup>®</sup> Curved Blades



### DYONICS ORBIT Curved Blades

Packaged 6 per box, sterile

Unique cutting window can be rotated up to 360° while in the joint

Allows for convex and concave blade orientation

Eliminates the need for removing and bending blades

Curve length 17 mm

	Power/EP-1° REF	Ultra Series REF†	Description	Latch Color
	7205320	7209029	4.5 mm ORBIT INCISOR <sup>®</sup>	Lime green
	7205321	7209030	4.5 mm ORBIT SYNOVATOR <sup>®</sup>	Forest green
	7205322	7209031	4.5 mm ORBIT RAZORCUT <sup>®</sup>	Pink
	7205323	7209028	4.5 mm ORBIT Full Radius	Yellow
	7210719	N/A	4.5 mm ORBIT Full Radius ELITE <sup>®</sup>	Maroon
	7205687	7209032	4.5 mm ORBIT INCISOR Plus	Violet

REF	Description
7205980	Cleaning brushes for curved blades, packaged 6 per box, not sterile

† Ultra Series blades are compatible with the DYONICS Power UltraLight Motor Drive only (REF 7205971 pg. D-2)

## Pre-Curved and Long Disposable Blades

### Pre-Curved Disposable Blades

All curve lengths 17 mm, packaged 6 per box, sterile

	Power/EP-1 REF	Ultra Series REF†	Description	Latch Color
	7211033	N/A	4.5 mm concave, curved, TURBO WHISKER® Blade	Navy blue
	7205332	7209026	4.5 mm concave, curved, INCISOR®	Lime green
	7205334	7209027	4.5 mm concave, SYNOVATOR®	Forest green
	7210980	N/A	3.5 mm concave, Full Radius	Beige
	7205349	7209025	4.5 mm concave, Full Radius, curve	Yellow
	7210521	N/A	4.5 mm concave, Full Radius with Incisor Inner, curve	Lilac
	7210988	N/A	4.5 mm concave, Full Radius ELITE®	Maroon
	72200036	N/A	4.5 mm INCISOR Plus ELITE	Slate
	7209075	N/A	4.5 mm concave, curved INCISOR® Plus	Violet

### Smith & Nephew DYONICS® POWER/EP-1® Long Shaver Blades and Burrs





18 cm working length, packaged 3 per box, sterile

Power/EP-1 REF	Description	Latch Color
7205335	4.5 mm curved, concave, Full Radius	Yellow
7205336	4.5 mm curved, convex, Full Radius	Yellow
7205341	4.5 mm straight, Full Radius	Yellow
7205343	4.5 mm straight, INCISOR®	Lime green
72200494	4.5 mm Curved INCISOR Plus ELITE	Slate
72200414	4.5 mm INCISOR Plus ELITE	Slate
6900678	5.5 mm Full Radius	Clear
72200080	4.0 mm Long ABRADER	Purple
72200082	5.5 mm Long ABRADER	Black

† Ultra Series blades are compatible with the DYONICS Power UltraLight Motor Drive only (REF 7205971 pg. D-2)

## Straight Disposable Blades

### Straight Disposable Blades Packaged 6 per box, sterile

	Power/EP-1 REF	Description	Latch Color
	7205316	4.5 mm TURBOWHISKER°	Navy blue
	7209303	4.5 mm Ultra Series TURBOWHISKER†	Navy blue
	7205312	3.5 mm INCISOR°	Medium grey
	7209489	4.0 mm INCISOR	Violet
	7205313	4.5 mm INCISOR	Lime green
	7205314	5.5 mm INCISOR	Mocha
	7205310	4.5 mm SYNOVATOR°	Forest green
	7205311	5.5 mm SYNOVATOR	Melon yellow
	7209077	4.5 mm SYNOVATOR outer with Full Radius Inner	Clear
	7205317	3.5 mm RAZORCUT°	Purple
	7205318	4.5 mm RAZORCUT	Pink
	7205319	5.5 mm RAZORCUT	Light pink
	7205308	3.5 mm Cutter	Red
	7205309	4.5 mm Cutter	Blue
	7205685	4.5 mm TURBOCUTTER° II	Sage Green

† Ultra Series blades are compatible with the DYONICS Power UltraLight Motor Drive only (REF 7205971 pg. D-2)

## Straight Disposable Blades

Straight Disposable Blades  
Packaged 6 per box, sterile

	Power/EP-1 REF	Description	Latch Color
	7205305	3.5 mm Full Radius	Beige
	7205306	4.5 mm Full Radius	Yellow
	7205307	5.5 mm Full Radius	Orange
	7209201	4.5 mm Full Radius Plus	Burnt orange
	7210751	3.5 mm Full Radius ELITE <sup>°</sup>	Beige
	7210499	4.5 mm Full Radius ELITE	Maroon
	7210909	5.5 mm Full Radius ELITE	Olive
	7205315	4.5 mm TURBOTRIMMER <sup>°</sup>	Powder blue
	7209202	4.5 mm TRIMMER	Brown
	7205345	4.5 mm INCISOR <sup>°</sup> PLUS	Violet
	7205459	5.5 mm INCISOR PLUS	Teal
	7220095	3.5 mm INCISOR PLUS ELITE	White
	7210976	4.5 mm INCISOR PLUS ELITE	Slate
	7220081	5.5 mm INCISOR PLUS ELITE	Lilac
	7209211	3.2 Chondroplasty Blade	

## Disposable Burrs

### Straight Disposable Burrs Packaged 6 per box, sterile

	Power/EP-1 REF (Wide Mouth)	Description	Latch Color	Maximum RPM
	7205324	4.0 mm Abrader	Aqua	8,000
	7205325*	5.5 mm Abrader	Black	8,000
	72200722	4.0 mm ELITE° Abrader	Aqua	10,000
	72200723	5.5 mm ELITE Abrader	Black	10,000
	7209706	4.0 mm Abrader, Left-handed	Clear	8,000
	7205326	4.0 mm Acromionizer	Mauve	8,000
	7205327*	5.5 mm Acromionizer	Brown	8,000
	7209512	4.0 mm Acromionizer with large suction window	Mauve	8,000
	72200724	4.0 mm ELITE Acromionizer	Mauve	10,000
	72200725	5.5 mm ELITE Acromionizer	Brown	10,000
	7205330	4.0 mm STONECUTTER®	Maroon	8,000
	7205331*	5.5 mm STONECUTTER	Olive	8,000
	72200728	4.0 mm ELITE STONECUTTER	Maroon	10,000
	72200729	5.5 mm ELITE STONECUTTER	Olive	10,000
	7209069	5.5 mm STONECUTTER, Left-handed	Clear	8,000
	7205328	4.0 mm NOTCHBLASTER®	Lilac	8,000
	7205329*	5.5 mm NOTCHBLASTER	Peach	8,000
	72200726	4.0 mm ELITE NOTCHBLASTER	Lilac	10,000
	72200727	5.5 mm ELITE NOTCHBLASTER	Peach	10,000
	7205668	4.0 mm ACROMIOBLASTER®	Sage green	8,000
	7205669*	5.5 mm ACROMIOBLASTER	Brick red	8,000
	72200730	4.0 mm ELITE ACROMIOBLASTER	Sage green	10,000
	72200731	5.5 mm ELITE ACROMIOBLASTER	Brick red	10,000
	7209287	5.5 mm ACROMIOBLASTER, Left-handed	Clear	8,000
	7209335	4.5 mm Shielded Burr	White	8,000
	6900988	4.0 mm Tapered Burr	Clear	8,000
	6900996	5.5 mm Tapered Burr	Clear	8,000
	7205727	4.5 mm HELICUT®	Slate	5,000
	7210460	4.5 mm HELICUT, Left-handed	Clear	5,000

\* Item shown

## Disposable Small Joint Blades/Burrs

Driven by the DYONICS II® Power Console, Smith & Nephew's latest platform for small joint arthroscopy includes updated hand controls, enhanced suction lever and port, completely redesigned blade hub attachment, and advanced performance.



REF	Blade Style
72201500*	DYONICS® POWERMINI® Small Joint Handpiece with blade multi-positioning, hand controls, suction lever, 10 ft. (3m) cord
72201503*	DYONICS® POWERMINI® Small Joint Handpiece with blade multi-positioning, no hand controls, suction lever, 10 ft. (3m) cord

	REF	Blade Style	Latch Color
	72201504	2.0 mm TURBOWHISKER®, 7 cm length	Yellow
	72201505	2.9 mm TURBOWHISKER, 7 cm length	Green
	72201506	2.9 mm Cutter, 7 cm length	Beige
	72201507	2.0 mm Full Radius, 7 cm length	Blue
	72201508	2.5 mm Full Radius, 7 cm length	Powder Blue
	72201509	2.9 mm Full Radius, 7 cm length	Red
	72201510	3.5 mm Full Radius, 7 cm length	Plum
	72201512	3.5 mm RAZORCUT®, 7 cm length	Pumpkin Orange
	72201514	3.5 mm INCISOR PLUS ELITE, 7 cm length	Sage Green
	72201518	2.9 mm Abrader, 7 cm length	Orange
	72201519	3.5 mm Abrader, 7 cm length	Aqua

## Disposable Mini Blades/Burrs



### Disposable Mini Blades/Burrs

Used with the Smith & Nephew DYONICS POWER<sup>®</sup>/EP-1<sup>®</sup> Mini Motor Drive (REF 7205357 pg. D-2).

Indications: Small joint arthroscopy. Packaged 6 per box, sterile

	Power/EP-1 REF	Description	Hub Color
	3606	2.9 mm Cutter, 7 cm length	Beige
	3407*	2.0 mm TURBOWHISKER <sup>®</sup> , 7 cm length	Yellow
	3418	2.9 mm TURBOWHISKER <sup>®</sup> , 7 cm length	Green
	3410*	2.0 mm Full Radius, 7 cm length	Blue
	3419	2.9 mm Full Radius, 7 cm length	Red
	6900927	3.5 mm Full Radius, BONECUTTER <sup>®</sup> , 7 cm length	Black
	3553*	2.9 mm Barrel Abrader, 7 cm length	Purple
	3530*	2.9 mm Abrader, 7 cm length	Orange
	3819	3.5 mm Abrader, 7 cm length	Aqua
	4224*	3.5 mm RAZORCUT <sup>®</sup> , 7 cm length	Pink
	4225*	3.5 mm INCISOR <sup>®</sup> , 7 cm length	Turquoise
	7209065	2.9 mm INCISOR, 7 cm length	Lime green
	7209085	2.0 mm INCISOR, 7 cm length	Navy blue

\*item shown

## Disposable Mini Blades — Long and Curved



Used with Smith & Nephew DYONICS® POWER/EP-1® Mini Motor Drives (REF 7205357 pg. D-2)





Curved blade allows greater access to anatomy

Blade Window Lock feature sets blade to lock at any position in window, as selected by the surgeon, protecting tissue from blade edges and controlling suction when not in cutting mode

Long and curved disposable mini blades are identical to the current mini blades, with one clear difference: they are one centimeter longer (8 cm)

### Disposable Mini Blades — Long and Curved

Indications: Small joint arthroscopy. Packaged 6 per box, sterile

	Power/EP-1 REF	Description	Hub Color
	7205474	3.5 mm Full Radius, 8 cm length	Plum
	7210488	3.5 mm Full Radius curved, 8 cm length concave	Brick red
	7205477	3.5 mm RAZORCUT®, 8 cm length	Pumpkin/orange
	7205558	3.5 mm curved, concave, INCISOR®, 8 cm length	Brick red
	7205476	3.5 mm INCISOR Plus, 8 cm length	Sage green

### Mini Irrigating Blades

Indications: Functional Endoscopic Sinus Surgery. Packaged 6 per box, sterile

Power/EP-1 REF	Description	Hub Color
7205688	3.5 mm Full Radius, irrigating, 9 cm length	Plum
7205692	3.5 mm Curved concave INCISOR, irrigating, 9 cm length	Brick red
7205693	3.5 mm INCISOR Plus, irrigating, 9 cm length	Sage green
7205694	FESS Blade Suction tubing	

## Short Sheath Components

### Short Sheath and Obturator

Packaged 6 per box, sterile

Power/EP-1 REF	Description	Latch Color
7205468	4.0 mm Short Sheath and Obturator handle	Maroon
7205470	5.5 mm Short Sheath and Obturator handle	Olive


### Short Sheath and Obturator Components

Reusable Obturator and Replacement O-rings: DYONICS® POWER/EP-1®.

Packaged 12 per box



Power/EP-1 REF	Description
7205458	Reusable O-rings, 4.0 mm
7205457	Reusable O-rings, 5.5 mm



7205456	Reusable obturator, 4.0 mm
7205455	Reusable obturator, 5.5 mm

