

N

Smith & Nephew 560P Camera
Control Units.....N-1

Smith & Nephew 450P/460P
Camera Control UnitsN-2

Smith & Nephew 470H Progressive
Scan Autoclavable Camera HeadN-3

Surgical Video

Product Highlights



560P HD Camera Control Unit
Providing exceptional versatility for
endoscopic procedures in all OR settings.

Smith & Nephew 560P Camera Control Units



Smith & Nephew 560P Camera Control Unit

Smith & Nephew's 560P Camera Control Unit (CCU) provides exceptional high definition image quality, broadcast quality color reproduction, and versatility for all endoscopic procedures in the OR setting. This new generation of high definition CCU will meet or exceed your facility's multi-specialty needs. Compatible with the latest in digital technology with full 1080 output in DVI, HD-SDI, and YPrPb with a powerful user menu for customization.

REF	Description
	Complete system includes: 560P Camera Control Unit; power cord; SDI/HDTV Cable, RG-59, 10', and country specific manual
72200571	560P Camera Control Unit, NTSC/PAL, U.S. power cord, 100-240 VAC, 50/60 Hz
72200572	560P Camera Control Unit, NTSC/PAL, Spain, 100-240 VAC, 50/60 Hz
72200573	560P Camera Control Unit, NTSC/PAL, France, 100-240 VAC, 50/60 Hz
72200574	560P Camera Control Unit, NTSC/PAL, Germany, 100-240 VAC, 50/60 Hz
72200575	560P Camera Control Unit, NTSC/PAL, Italy, 100-240 VAC, 50/60 Hz
72200576	560P Camera Control Unit, NTSC/PAL, Sweden, 100-240 VAC, 50/60 Hz
72200577	560P Camera Control Unit, NTSC/PAL, Dutch, 100-240 VAC, 50/60 Hz
72200578	560P Camera Control Unit, NTSC/PAL, U.K., 100-240 VAC, 50/60 Hz
72200579	560P Camera Control Unit, NTSC/PAL, International, 100-240 VAC, 50/60 Hz
72201380	560P Camera Control Unit, NTSC/PAL, Korea, 100-240 VAC, 50/60 Hz
72201381	560P Camera Control Unit, NTSC/PAL, Portugal, 100-240 VAC, 50/60 Hz
72201382	560P Camera Control Unit, NTSC/PAL, Norwegian, 100-240 VAC, 50/60 Hz
72201383	560P Camera Control Unit, NTSC/PAL, Danish, 100-240 VAC, 50/60 Hz



Smith & Nephew 560H Non-autoclavable Camera Head

Smith & Nephew's 560H broadcast high definition camera head is compatible with a complete line of Smith & Nephew direct-view endoscope, HD videoarthroscopes, monitors, and 560P HD Camera Control Unit. Provides for 3 buttons (6 function) control of CCU, peripheral devices, printers and Smith & Nephew Image Management Systems.

REF	Description
72200561	560H Camera Head, HD, NTSC/PAL

Smith & Nephew 450P/460P Camera Control Units



Smith & Nephew 450P/460P Camera Control Units

Smith & Nephew's 450P and 460P Camera Control Units (CCUs) provide exceptional versatility for endoscopic procedures in all OR settings. This new generation of CCUs will meet your facility's multi-specialty needs by providing procedure-specific image optimization as well as support for both autoclavable and non-autoclavable solutions. Compatibility with progressive scan camera heads and support for future software upgrades help ensure your investment remains secure and up-to-date.

In addition, the 460P CCU includes the latest in digital technology for state-of-the-art image fidelity and communications, including a complete set of digital video outputs (DVI, SDI, and DV) and digital zoom capability.

REF	Description
	Complete system includes: Camera Control Unit; power cord; 1 6-pin-to-6-pin, 15 ft IEEE-1394 interconnected cable (REF 8109734) (460P only); 1 RGB cable (REF 130-0068); 1 S-video Cable (REF 130-0053); country specific manual
72200085	450P Camera Control Unit, NTSC, U.S. power cord, 120 VAC
72200086	460P Camera Control Unit, NTSC, U.S. power cord, 120 VAC
72200103	450P Camera Control Unit, PAL, U.K. power cord, 230 VAC
72200104	460P Camera Control Unit, PAL, U.K. power cord, 230 VAC
72200105	450P Camera Control Unit, PAL, International power cord, 230 VAC
72200106	460P Camera Control Unit, PAL, International power cord, 230 VAC

Smith & Nephew 460H Non-autoclavable Camera Head



Smith & Nephew 460H camera heads are compatible with a complete line of Smith & Nephew direct-view endoscopes, VideoEndoscopes (C-Mount), monitors, and camera control units. Three buttons allow remote control of compatible peripheral devices like the CCUs, printers, and Smith & Nephew Image Management Systems.

The Smith & Nephew 460H digital 3 CCD camera head is engineered to be a plug and play imaging solution and to maximize hospital efficiency. From the very creation of the image to the surgical environment, new technology in the 460H camera head successfully delivers all the aspects of a high quality surgical image for every endoscopic procedure.

REF	Description
72200090	460H Digital Camera Head, NTSC
72200091	460H Digital Camera Head, PAL

Smith & Nephew 470H Progressive Scan Autoclavable Camera Head



Smith & Nephew 470H Progressive Scan Autoclavable Camera Head

Smith & Nephew 470H camera heads are compatible with a complete line of Smith & Nephew direct-view endoscopes, VideoArthroscopes, monitors, and camera control units. Three buttons allow remote control of compatible peripheral devices like the CCUs, printers, and Smith & Nephew Image Management Systems.

The Smith & Nephew 470H Autoclavable 3 CCD camera head is engineered to be a plug and play imaging solution and to maximize hospital efficiency. The 470H head delivers exceptional image quality in an extremely durable and ergonomic package built to withstand the intense environment of the autoclave.

REF	Description
72200092	470H Autoclavable Camera Head, NTSC
72200093	470H Autoclavable Camera Head, PAL

Replacement Parts

81010558	Replacement kit (soak cap, radial and crush seal O-rings, and lanyard)
8109222	O Rings
8109223	O Rings

