

O

Video Image Capture Devices.....	O-1
Video Monitors.....	O-2
Video Printers and Recorders.....	O-3
Video Printer Paper/Printer Accessories....	O-4
Video Couplers and Adaptors.....	O-5
Video and Accessory Cables.....	O-7
Video Carts and Accessories.....	O-8
Made to Order.....	O-9

Video Accessories



660HD Image Management System

Image management increases OR efficiency by easily capturing and utilizing surgical images in true high definition (HD) format.

Video Image Capture Devices

660HD Image Management System

The Smith & Nephew 660HD Image Management System is an easy-to-use digital capture system that enables the OR staff to capture, edit, print, and store intraoperative surgical images and videos in true 1920 x 1080 high definition (HD) format.



CONDOR-Enabled 660HD Image Management System

REF	Description
72200812	660HD Image Management System, English on-screen display and Operations/Service Manual, U.S. power cord
72200814	660HD Image Management System, multi-language on-screen display and Operations/Service Manual, continental power cord
72200816	660HD Image Management System, English on-screen display and Operations/Service Manual, U.K. power cord

Accessories

7210344	Touch panel, LCD, 15" (domestic)
7210345	Touch panel, LCD, 19" (domestic)
7209339	Single-pedal footswitch
7209340	Dual-pedal footswitch
7210534	DVD-R media
7209338	650 MB or 700 MB CD-R media
72201464	Printer HP7460
7210533	Premium photo paper, high-gloss, 50 sheets
72201745	Premium photo paper, soft-gloss, 100 sheets
72200204	HP 99 Photo Inkjet Print Cartridge with Vivera Ink
72200205	HP 97 Tri-color Inkjet Print Cartridge, high volume

Video Monitors



Smith & Nephew flat panel displays are optimized for use with the full range of Smith & Nephew video cameras and are part of the recommended Smith & Nephew visualization system configurations. These flat panel displays offer exceptionally clean, crisp, high-resolution images, ensuring superior reproduction of anatomical detail. Full HD widescreen displays deliver the full resolution of the Smith & Nephew 560 HD camera, 1920x1080i. HD compatible displays can display images from standard definition video sources but have lower than HD resolution.

HD Widescreen, 16:9, Flat Panel Displays

REF	Description
72201421	24" HD Widescreen HD 1080i LCD
72201422	24" HD Widescreen HD 1080i LCD, UK
72201423	24" HD Widescreen HD 1080i LCD, International

HD Compatible Flat Panel Displays

REF	Description
72201434	19" HD Compatible 5:4 square LCD
72201435	19" HD Compatible 5:4 square LCD, UK
72201436	19" HD Compatible 5:4 square LCD, International
72201085	32" HD Compatible Widescreen HD LCD
72201086	32" HD Compatible Widescreen HD LCD, UK
72201087	32" HD Compatible Widescreen HD LCD, International

Desktop Monitor Stand

REF	Description
72200019	Desktop Monitor Stand for LCDs below 25 inches, Sony® SU-560

Video Printers and Recorders

Video Printers and Recorders

REF	Description
72200992	UP-55MD A5 Color Video Printer by Sony®, NTSC/PAL, RGB, S-VHS, composite video, compatible with Smith & Nephew digital camera head control buttons, medical grade, 8MB memory, 100–240 VAC
72201464	Plain paper printer. Up to 6 pages per minute in normal color mode. 30 pages per minute in draft color mode. Optional high-capacity ink cartridges so fewer cartridge changes are needed. 16 mb memory. Up to 4800-optimized dpi for professional photo quality. Advanced control panel with ink status display, and print quality selector.
7210817	Medical Grade DVD Video Recorder by Panasonic®. Real-time recording to digital video-in-video format, DVD-R, DVD-RAM compatible, controllable using a single-pedal footswitch or CCU accessory inputs, can create DVDs that play on commercial DVD players or PCs.

Video Printer Paper/Printer Accessories

Video Printer Paper

REF	Description
72200993	Color printing pack for Sony® UP-55MD printer (REF 72200992), 1 pack includes 2 bags of printer paper (100 each) and 2 color ribbon refills.
7207802	Sony® part # UPC-510, color printing pack for Sony® UP-51MD/MDU series printers (Smith & Nephew REF 7207801); 1 pack includes 2 bags of print paper (100 each) and ink ribbon; 200 prints per pack
7204983	Sony® part # UPC-5510, color printing pack for Sony® UP-5600MD series printers (Smith & Nephew REF 7204973); 1 pack includes 2 bags of print paper (100 each) and ink ribbon; 200 prints per pack
7206058	Sony® part # UPC-5510C, color printing pack for Sony® UP-5600MD series printers (Smith & Nephew REF 7204973); 1 pack includes 2 bags of print paper (100 each) and 1 disposable ribbon-loaded cartridge; 200 prints per pack
7206097	Sony® part # UPC-2010, color printing pack for Sony® UP-2900MD series printers (Smith & Nephew REF 7206094); 1 pack includes print paper and ribbon; 200 prints per pack
4249	Video Print Papers/Ink Cartridge for Toshiba printer (Smith & Nephew REF 7204253); 100 prints
3698	Sony® part # UPC-5010A, Video Print Papers/Ink Cartridge for Sony® UP-5000 series printers (Smith & Nephew REF 3758, 4238, and 4259); 100 prints
4217	Sony® part # UPC-3010 Video Print Papers/Ink Cartridge for Sony® 3000 series printers (Smith & Nephew REF 4150); 80 prints
7205206	Sony® part # UPC-1010 Video Print Papers/Ink Cartridge for Sony® 1800 series printers (Smith & Nephew REF 7204998); 100 prints
7204994	Replacement ribbon holder for Sony® UP-5600 series printers (Smith & Nephew REF 7204973)

Smith & Nephew, Inc. is a recognized Sony® distributor. All Sony® products not shown in this catalog are available upon request.

Printer Accessories

7210533	HP Premium Plus High Gloss Photo Paper, 8.5" x 11", 50 sheets per package, for use with REF 72200035
72200204	HP High Volume Ink Cartridge, black 21 ml, for use with REF 72200035
72200205	HP High Volume Ink Cartridge, color 14 ml, for use with REF 72200035
72101745	HP Premium Plus Soft-gloss Photo Paper, 8.5" x 11", 50 sheets per package, for use with REF 72200426

Video Couplers and Adaptors

Non-Autoclavable Endocouplers

Superior clarity and brightness

Optimal picture quality

Ease-of-use: convenient one-finger scope attachment

Long-term reliability

Universal application: couplers fit rigid eyepiece endoscopes to V- or C-Mount camera heads of virtually any make or model

Compatible with EtO, Steris®, and Sterrad® sterilization modalities



REF

Description

72201635

Coupler, High Definition, Non-autoclavable, 19.5 mm, Black ring

72201730

Coupler, High Definition, Non-autoclavable, 26.0mm, Silver ring



REF

Description

7204613

C-Mount, 25 mm Coupler (black focus ring)

7204823*

C-Mount, 30 mm Coupler (gold focus ring)

7204614

C-Mount, 35 mm Coupler (blue focus ring)



7205461

C-Mount Zoom Coupler, 25-45 mm, for use with Smith & Nephew Digital 3-Chip Camera Models D3, ED-3, and HD 460H



4533

C-Mount, 35 mm, 90° Urology Coupler, free-floating (lockable), rotating head



4464

C-Mount, 25 mm Beamsplitter, free-floating (lockable), rotating head

Replacement Parts

REF

Description

8001701

Replacement O-ring for Videoscopes and C-Mount coupler

8012760

Replacement O-ring for V-Mount coupler

* Item shown

Video Couplers and Adaptors

Autoclavable Couplers

Superior clarity and brightness

Thread compatible with Smith & Nephew autoclavable camera offering

Unparalleled brightness

Sharp images throughout the entire field-of-view

Sapphire lenses for improved scratch resistance over conventional glass

Epoxy-free, proprietary seal design guarantees leak-free autoclave performance

Compatible with steam and chemical (EtO, Steris®, and Sterrad®) sterilization modalities

High Definition Endocouplers – 560H head only



REF	Description
72200315	19.5 mm High Definition, C-mount, Black ring
72200422	26.5 mm High Definition, C-mount, Silver ring

Standard Definition Endocouplers



REF	Description
7210192	18 mm Smith & Nephew mount for use with 370 and 470H Camera Heads only

Video and Accessory Cables

Video and Accessory Cables

REF	Description
7206083	Accessory Control Cable, mini-phono to mini-phono, for use with 325Z, ED-3, and HD900 Video Cameras for remote connection with Video Printers (Smith & Nephew REF 7207812 and 7206094), Video Cassette Recorder, and Smith & Nephew 600-series Image Management System
4154	Accessory Control Cable, RCA to mini-phono, 6 ft., for use with D3, D1, and MiniCam Video Cameras for remote connection with Video Printers (Smith & Nephew REF 7207801, 7207812, and 7206094), Video Cassette Recorder (REF 7205010), and Smith & Nephew 600-series Image Management System
7211032	Cable Remote Adaptor to enable control of the DVD Video Recorder (7210817) with a single-pedal footswitch or CCU accessory outputs.
4175	RGB Cable, 6 ft (1.8 m)
4174	RGB Cable, 25 ft (7.6 m)
7209718	RGB Cable, 50 ft (15.2 m)
7209719	RGB Cable, 75 ft (22.8 m)
7210114	RGB Cable with single connector for flat panel displays, 6 ft (1.8 m).
7204585	DB9 to DB25 RS232 Data Cable, for use with Smith & Nephew Digital Camera Control Unit and Video Print Recorders (REF 7204253, 4259, 4238, 4251)
72201420	SDI/HDTV Cable, RG-59, BNC, 10'
72201438	SDI/HDTV Cable, RG-59, BNC, 25'
7209709	High-Grade S-VHS (S-video) Cable, 1.5 ft (46 cm)
7204580	High-Grade S-VHS (S-video) Cable, 6 ft (1.8 m)
7204582	High-Grade S-VHS (S-video) Cable, 15 ft (4.5 m)
7204583	High-Grade S-VHS (S-video) Cable, 25 ft (7.6 m)
7204584	High-Grade S-VHS (S-video) Cable, 50 ft (15.2 m)
7209710	High-Grade S-VHS (S-video) Cable, 75 ft (22.8 m)
7205002	DB9-RCA Remote Cable
7205196	Talley Cable for Sony® SVO-2000
7204586	RCA to Mini Cable
7209761	S-USB Cable, 50 ft (15.2 m)
7209762	S-USB Cable, 75 ft (22.8 m)
7209775	S-USB Cable, 3 ft (.92 m)
7209836	S-USB Cable, 10 ft (3 m)
7209848	Audio Hook-Up Cable, 1/8" phono to 1/8" phono, 50 ft (15.2 m)
7209849	Audio Hook-Up Cable, 1/8" phono to 1/8" phono, 75 ft (22.8 m)
7209850	Audio Hook-Up Cable, 1/8" stereo phono to 1/8" stereo phono, 50 ft (15.2 m)
7209851	Audio Hook-Up Cable, 1/8" stereo phono to 1/8" stereo phono, 75 ft (22.8 m)
7209852	Audio Hook-Up Cable, 1/8" stereo phono male to 1/8" stereo phono female, 50 ft (15.2 m)
7209853	Audio Hook-Up Cable, 1/8" stereo phono male to 1/8" stereo phono female, 75 ft (22.8 m)
7209856	Audio Hook-Up Cable, 1/8" phono to 1/8" phono, 1.5 ft (46 cm)

Video Carts and Accessories



Endoscopy Base Carts

REF	Description
72202161	Cross Specialty Cart with LCD Video Boom
72202162	Cross Specialty Cart with LCD Video Boom and Transformer
72202163	Cross Specialty Cart with LCD Post-Mount
72202164	Cross Specialty Cart with LCD Post-Mount and Transformer

Made to Order

Made to Order products have longer production times. Please contact your sales representative or customer service for more information.

Urology Coupler



REF	Description
7209042	C-Mount, 25 mm, 90° Urology Coupler, free-floating (lockable), rotating head, for 4.0 mm and 5.0 mm endoscopes

