

R

Powered Instruments

Smith & Nephew DYONICS® POWER Instrumentation (battery-operated)	R-1
Accessories	R-2
Smith & Nephew DYONICS® POWER Instrumentation (corded)	R-4
Accessories	R-6

Smith & Nephew DYONICS[◇] POWER Instrumentation (battery-operated)

Battery Operated Power Equipment

REF	Description
2851	450 Drill, with battery, 1/4" Jacobs Chuck
	
2855	Sagittal Saw, with battery
	
2856	800 Wire Driver, with battery
	

Accessories

Charger, Power Pack, and Chuck Key







	REF	Description
	2863	Pac II Battery, 6.0 VDC, box of 2
	3082	MPC IV Charger, 110 VAC*
	3438	Jacobs Chuck Key

* Batteries not included with charger

Accessories

Sterile, Disposable Saw Blades†

For use with Sagittal Saw, REF 2855

	REF	Description
	3701	Saw Blade – Standard, thickness: 0.65 mm, length: 63 mm, box of 6
	3702	Saw Blade – Standard, thickness: 0.89 mm, length: 63 mm, box of 6
	3703	Saw Blade – Short, thickness: 0.65 mm, length: 29 mm, box of 6
	3704	Saw Blade – Short, Narrow, thickness: 0.65 mm, length: 29 mm, box of 6
	3705	Saw Blade – Narrow, thickness: 0.65 mm, length: 63 mm, box of 6
	3706	Saw Blade – Long, special thickness: 0.89 mm, length: 68 mm, box of 6

† All saw blades have graduated markings every 5 mm

Smith & Nephew DYONICS[®] POWER Instrumentation (corded)

Corded Power Equipment

For use with Smith & Nephew DYONICS[®] POWER Control Unit (see pg. D-2). The drill's design and choice of modular accessories provide smooth switching of components during arthroscopic surgery. For the sagittal saw, you may choose from multiple blades with varying sizes.

	REF	Description
	7205785	Pistol Grip Drill (requires REF 7205788)
	7205786	Inline Sagittal Saw (requires REF 7205788)
	7205791	Sagittal Saw Module (used with REF 7205785)
	7205787	Sterilization Tray (holds all components including drill and saw)

Smith & Nephew DYONICS[◇] POWER Instrumentation (corded)








Instrumentation

	REF	Description
	7209391	Smith & Nephew DYONICS POWER high-speed drill (model 6100)
	7209394	Smith & Nephew DYONICS POWER oscillating saw, MicroAire [®] interface (model 6710)
	7209392	Smith & Nephew DYONICS POWER reciprocating saw, MicroAire interface (model 6420)
	7209395	Smith & Nephew DYONICS POWER heavy-duty oscillating saw, MicroAire interface (model 6740)
	8103998	Smith & Nephew DYONICS POWER locking wrench
	7209393	Smith & Nephew DYONICS POWER Synthes chuck

Accessories

Wire Driver, Pin Driver, Chucks, and Cable



For use with Smith & Nephew DYONICS® POWER Control Unit.

	REF	Description
	7205789	Automatic Wire Driver (accepts wires up to 1.6 mm)
	7205790	Automatic Pin Collet (accepts pins up to 3.2 mm)
	7205792	Jacobs Chuck, 5/32"
	8104907	Chuck Key
	7205793	Jacobs Chuck, Keyless, 1/8"
	3438	Chuck Key
	7205794	Jacobs Chuck, 1/4"
	7205795	Jacobs Chuck, Keyless, 1/4"
	7205788	Power Cable (from drill and/or sagittal saw to control unit)

Accessories


Disposable Saw Blades for Smith & Nephew DYONICS® POWER Instrumentation

For use with Sagittal Saws REF 7205786 or 7205791. Disposable, packaged 6 per box, sterile.

	REF	Description
	7209396	DYONICS sagittal saw blade, 5.8 mm x 12 mm x 0.4 mm, with straight fine teeth, no etchings, depth stop
	7209397	DYONICS sagittal saw blade, 5.8 mm x 18 mm x 0.4 mm, with straight fine teeth with markings, depth stop





Sagittal Saw Blades (WxLxT)

For use with Sagittal Saws REF 7205786 or 7205791. Disposable, packaged 6 per box, sterile.

	7205809	Saw Blade, 9.4 mm x 15 mm x 0.4 mm, with 15 mm depth stop, Titanium Nitride finish, No etching
-----------------------------------------------------------------------------------	---------	------------------------------------------------------------------------------------------------

Sagittal Saw Blades with Depth Stops









For use with Sagittal Saws REF 7205786 or 7205791. Disposable, packaged 6 per box, sterile.

	7205807	Saw Blade, 9.4 mm x 10 mm x 0.4 mm, with 10 mm depth stop, Titanium Nitride finish, No etching
	7205808	Saw Blade, 9.4 mm x 12 mm x 0.4 mm, with 12 mm depth stop, Titanium Nitride finish, No etching
	7205898	Saw Blade, 9.4 mm x 10 mm x 0.4 mm, with 10 mm depth stop, without Titanium Nitride finish, No etching
	7205899	Saw Blade, 9.4 mm x 12 mm x 0.4 mm, with 12 mm depth stop, without Titanium Nitride finish, No etching
	7205900	Saw Blade, 9.4 mm x 15 mm x 0.4 mm, with 15 mm depth stop, without Titanium Nitride finish, No etching

Accessories

Calibrated Sagittal Saw Blades (WxLxT)

For use with Sagittal Saws REF 7205786 or 7205791. Disposable, packaged 6 per box, sterile.

	REF	Description
	7205802	Saw Blade, 19.1 mm x 25.4 mm x 0.4 mm, with markings
	7205803	Saw Blade, 5.8 mm x 25.4 mm x 0.4 mm, with markings
	7205804	Saw Blade, fine teeth, 9.4 mm x 25.5 mm x 0.4 mm, with markings
	7205805	Saw Blade, 9.4 mm x 25.5 mm x 0.4 mm, with markings
	7205806	Saw Blade, 4.5 mm x 25.5 mm x 0.4 mm, with markings
	7205895	Saw Blade, 19.1 mm x 25.4 mm x 0.4 mm, with markings and Titanium Nitride finish
	7205896	Saw Blade, 5.8 mm x 25.4 mm x 0.4 mm, with markings and Titanium Nitride finish
	7205897	Saw Blade, fine teeth, 9.4 mm x 25.5 mm x 0.4 mm, with markings and Titanium Nitride finish